



Key features:

- High frequency resolution for more precise tuning
- Automatic directionality for better hearing in noise
- Wide range of designs



The myRemote app enables easy control of default settings such as program selection and volume via smartphone.

iPhone is a registered trademark of Apple Inc. in the U.S. and other countries. Android is a registered trademark of Google Inc. The Bluetooth® word mark and logos are the property of Bluetooth SIG Inc. Any use of such marks by WSAUD A/S is governed by a license agreement.

Other trademarks and trade names remain the property of their respective owners.

08/2021



## Better hearing for everyone

P12 hearing aids are proven to improve speech understanding



### P12

Proven 12 compression channel technology



When you are looking for a professionally fitted yet affordable solution to improve your hearing, proven technology is a must. Our latest P12 hearing aids feature 50% more compression channels and enhanced directionality to improve your ability to hear what people are saying, even in noisy situations. With a wide range of designs, improved speech understanding is now more accessible than ever.

# Features for better hearing

Conversation made easy with Directional iLock



In places with strong background noise (e.g. traffic or family gatherings), communicating can require a lot of concentration and energy. Directional iLock helps to relieve this by focusing the sound on the person in front of you, reducing unwanted background chatter. The result is a more clear listening experience.

Phone calls made easy with XPhone



In noisy offices with phones ringing and co-workers are chatting, hearing during a business call can be challenging for anyone. The XPhone program simultaneously sends phone signals to both hearing aids when the phone is held to one ear, because we want you to hear the conversations clearly, even in a noisy office.

Increased spatial awareness



P12 hearing aids work to mimic the human ear and natural brain function to ensure you have good awareness of your surroundings. This is achieved by a combination of transferring signals between the left and right hearing aids, while applying a directional effect that helps make your environment sound natural.

# A wide range of models for everyone

**RIC**  
Flexibility and practicality in a compact package



Colors: Beige, Gray, Granite



**IF**  
Instant Fit and nearly invisible

Colors: Shell (Red, Blue, Mocha), Faceplate

**Custom**  
Perfectly fitted for ultimate comfort



Colors: Beige, Brown, Tan

# BTE models – the classic design

**M**  
Compact yet powerful



Colors: Beige, Gray, Granite

**P**  
Flexible and intuitive



Colors: Beige, Gray, Dark Brown

**XP**  
High power performance



Colors: Beige, Gray, Dark Brown

# Accessories for better hearing

**Smart Connect**  
Remote control and audio streaming device in one



**Smart Transmitter**  
Connects television, stereo systems and other devices without Bluetooth function with Smart Connect



**Smart Remote**  
Easily control your hearing aids remotely

